THE UNIVERSITY OF QUEENSLAND
DEPARTMENT OF PHYSICS
PHYS2041
CAVENDISH EXPERIMENT

Introduction

In this laboratory you will repeat the experiment in which Newton's gravitation constant, "G", was first measured. This constant is usually called "Big G", and the local acceleration of gravity is called "Little g". This experiment has also been called "weighing the earth" (why?). Although it falls outside the realms of quantum physics, it is included in this laboratory as attempts are still ongoing to improve the accuracy of the measurement. Currently G is only known to around 1 part in 10,000 and it is one of the least accurately known physical constants, despite being one of the first measured.

The apparatus consists of a dumbbell suspended on a very weak torsion fibre in the shape of a *tape* 10 microns thick and 150 microns wide. When large masses are brought near the small masses at the ends of the dumbbell in the horizontal plane, the dumbbell masses will rotate toward them until the torsion fibre just balances the gravitation force. In order to increase the sensitivity of the experiment, the torsion fibre must be very weak in order to allow an appreciable rotation of the dumbbell. However, this means that the period of the resultant oscillation will be very long. In our particular case, the oscillation period is about 10 minutes. Therefore it is very important that the apparatus is not bumped or touched during the course of the experiment. Any extraneous vibration takes a long time to die out. Also, any large (>70kg) masses should be kept away from the balance for the duration of the experiment. *This is a very delicate measurement!*

Experiment

At time-zero, the large masses are *very carefully* moved from one extreme orientation to the other so that they are *just touching* the glass cover. The position of the light spot on the scale should be read, recorded and plotted at 30 second intervals for about five oscillations of the dumbbell (ie for about 50 minutes). This will allow you to determine the torsion constant of the fibre as well as the new equilibrium position of the dumbbell. Next, move the large masses to their original extreme positions and continue the measurements for another five periods.

You should plot the data *as you go* in order to be able to see any problems which may occur, and in order to help you decide how you will extract the required information.

It is important to consider the effect of *both* large masses on *each* of the small masses. Also, you should take into account the damping of the motion when determining the torsion constant.

There are many small corrections of the order of a percent or so, but you can derive an exact solution by first finding the gravitational potential and then differentiating to obtain the torque.

The physical data on the apparatus supplied by the manufacturer are


| Mass (g) | Radius (mm) |
|----------|-------------|
| 15 | 6.9 |
| 1500 | 32 |

The dumbbell radius is 50 mm and the distance between the centre of the large ball and the centre of the chamber when the large ball is touching the glass is 46.5 mm.

Reference

Physics for Scientists and Engineers, Tipler PA, Freeman & Co.

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$$f = \frac{S \cdot d \cdot b^2}{2 m_1 \cdot t^2 \cdot L} \quad (II)$$

The device can be used to demonstrate the gravitational attractive force between masses and to determine the gravitational constant f .

Measurement principle (Fig.1)

The static equilibrium (position I) of a torsion pendulum with a period of approx. 10 min. is disturbed by a change in position of the outer masses m_1 , which affect the dumbbell-shaped pendulum body. The oscillations become damped and the pendulum takes up a new equilibrium position (position 2). The angle between the two positions is a measure of the active gravitational force.

The oscillation of the pendulum, which is equipped with a concave mirror, is indicated by a light pointer. This is possible using

- visible light directly on a mm-scale or
- infrared light along with the IR position detector (332 11). The latter enables the measurement values to be plotted or evaluated by computer.

The gravitational constant f can be obtained from the oscillation curve with respect to time, the mass m_1 and the geometry of the arrangement using either the end deflection method or (in a quicker process) the acceleration method.

In the *end deflection method* the torsion pendulum period T and the distance S separating the light pointer positions are evaluated for the two equilibrium positions:

$$f = \frac{\pi^2 \cdot b^2 \cdot d \cdot S}{m_1 \cdot T^2 \cdot L} \quad (I)$$

In the *acceleration method*, the acceleration of the torsion pendulum $a = 2s/t^2$ is evaluated using the masses m_1 after the equilibrium position has been disturbed:

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1 Safety notes

Protect the sensitive bronze band of the torsion pendulum from uncontrolled mechanical loading:

- Do not unscrew the oscillating system's locking screws ⑦ - see Fig. 2 - before the device has been correctly assembled and brought into position.
- Always lock the oscillating system when the device is not in use. In particular, make sure it is locked during transport and assembly.

The knurled screw (3.1) - see Fig. 2 - for fixing the torsion head (preadjusted on delivery) is only to be slightly loosened if a fine adjustment of the zero point proves necessary on putting the correctly-assembled torsion balance into operation. The grub screw (3.2) - see Fig. 2 - for fixing the pendulum holder is only to be loosened when replacing the torsion band as described in section 4.

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2 Description, technical data, scope of supply (see Fig. 2)


Fig.2
Gravitations-Drehwaage.
Gravitation torsion balance

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- ① Metal housing (15 cm dia.), with metal ring and glass covers; two-part slide (1.1) to prevent interference from convection currents between the housing and protective tube ②
- ② Protective tube (25 cm long) for torsion band ④
- ③ Torsion head with pendulum holder, can be rotated for adjustment of the pendulum equilibrium position when screw has been loosened (3.1); pendulum holder fixed with grub screw (3.2)
Important!
Do not loosen the knurled screw (3.1) unless it is essential to carry out a zero-point adjustment in accordance with section 3.5. Only loosen the grub screw (3.2) when the torsion band is being replaced as described in section 4.
- ④ Torsion band made from bronze, 26 cm long
Spare part no.: 683 21
- ⑤ Dumbbell-shaped pendulum body, consisting of 2 lead balls (5.1) on a metal rod

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- ⑥ Concave mirror for light pointer indication of pendulum motion
Focal length f approx. 30 cm
- ⑦ Screws for locking the pendulum body ⑤ using pair of springs (7.1)
- ⑧ Pair of large lead balls
Spare part no. for 1 ball: 683 22
- ⑨ Ball carrier, can be rotated around stand rod ⑪, for correct experimental arrangement of the large lead balls ⑧
- ⑩ Supporting ring with fixing screw (10.1) for ball carrier
- ⑪ Stand rod (9 cm x 1.2 cm dia.) for assembly of the device in stand material

Included in scope of supply:

1 m self-adhesive scale tape with cm and mm divisions

Important experimental data (see Fig. 1):

Housing depth: 30 mm

Torsion pendulum

Period: approx. 10 min

Diameter / mass m_2 of a lead ball (5.1):
15 mm / 20 g

Distance d between the center of a ball and the axis of rotation: 50 mm

Diameter / mass m_1 of a large lead ball ⑧:
64 mm / 1.5 kg ± 5 g

Distance b between the center of the large ball (when contact is made with housing) and the small ball (in the equilibrium position): 47 mm

3 Operation

Important!

Satisfactory experiment results are only possible when the torsion pendulum has been properly adjusted and the torsion oscillations produced by attraction between the masses are not affected by unwanted pendulum movements. The pendulum is very sensitive to any disturbance of the experiment setup: make sure that the experiment setup is absolutely stable, e.g. by attaching it to a solid wall or placing it on a sturdy bench or table.

Temperature variations cause convection in the housing of the torsion balance, which in turn cause undesired motions of the torsion pendulum (see Fig. 6).

For this reason, select an experiment site which does not stand in direct sunlight or drafts.

3.1 Additional equipment required

3.1.1. Oscillation indication using (visible) light marks on a mm-scale

| | |
|--|--------|
| 1 Lamp housing | 450 60 |
| 1 Lamp, 6 V, 30 W | 450 51 |
| 1 Aspherical condenser (with 1-mm slit screen) | 460 20 |
| 1 A.C. voltage source, 6 V, 30 W e.g. | 562 73 |

| | |
|----------------------|--------|
| 1 Steel tape measure | 311 77 |
| 1 Stop-clock e.g. | 313 05 |

| | |
|--------------------|--------|
| 1 Suspended magnet | 510 44 |
|--------------------|--------|

Wall assembly (see Fig. 3)

| | |
|--|--------|
| 1 Large stand base | 300 31 |
| 1 Pair of levelling screws | 300 06 |
| 1 Leybold multiclamp | 301 01 |
| 1 Rotatable clamp | 301 03 |
| 1 Stand rod, 47 cm | 300 42 |
| Assembly aids and accessories: Electric drill, masonry drill bit, wall plugs (6 mm) | |

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3.1.2 Recorder plotting or computer-assisted recording of the oscillation using infrared light; rail assembly on a bench (see Fig. 4)

| | |
|---|--------------|
| 1 Infrared position detector | 332 11 |
| 1 A.C. voltage source, 12 V e.g. | 562 73 |
| 1 TY-recorder e.g. | 575 701 |
| oder | |
| MS DOS computer via serial interface | |
| RS 232 cable | e.g. 530 008 |
| Program with experiment-specific default settings e.g. contained on the demo disk for "Universal Data Acquisition", included in the scope of supply of the IR position detector (332 11); see measuring example Fig. 4.1. | |
| oder | |
| MS DOS computer with interface | |
| CASSYpack-E | 524 007 |
| Program "Measuring and Evaluating" | 524 112 |
| 1 Steel tape measure | 311 77 |
| 1 Suspended magnet | 510 44 |
| 1 Optical bench with normal profile, 1 m | 460 32 |
| 2 Optical riders e.g. | 460 351 |
| 1 Stand rod, 25 cm e.g. | 300 41 |

Note:
 Further information on the use of the infrared position detector can be found in the corresponding Instruction Sheet.

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3.2 Mounting the stand base on the wall (required for wall mounting as shown in Fig. 3 before using the setup for the first time.)

Important:

Wall mounting of the stand material is not suitable for the experiment setup with the IR position detector 332 11 (see section 3.1.2; Fig. 4).

Correctly mount the stand base on its apex hole as shown in Fig. 5.1 - 3.3 and align parallel to the wall using levelling screws; use solid base material if the levelling screws press into the wall.


Fig. 5.1-5.3

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3.3 Preparing for the experiment

Assemble the arrangement - initially without the large lead balls - on the wall as shown in Fig. 3 (or similarly on the optical bench on a sturdy lab bench).

Align the lamp housing so that the reflected beam is projected unobstructed onto the scale, which should be positioned at least 5 m distant.

Align the filament so that it is vertical before attaching the condenser with screen holder and slit (vertically!) to the lamp housing.

Move the insert in the lamp housing so that a focused image of the filament is projected onto the torsion balance mirror. To do this, hold a white sheet of paper, on which the filament forms an image, directly in front of the balance.

Move the lamp housing so that the slit is focused onto the scale.

If necessary, adjust the zero-point as described in section 3.5.

Attach the lead balls and bring the arrangement to an extreme position. In carrying out this step make absolutely sure that the housing does not come into contact with either your finger or a lead ball.

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3.4 Experiment procedure; measuring example

Important!

After the experiment has been set up as described in section 3.3, it must be left undisturbed for at least two hours so that the pendulum can settle in its equilibrium state.

Turn the ball carrier ⑨ without touching the housing.

To shorten this settling time, make use of the diamagnetism of the lead: each time one of the pendulum lead balls approaches the glass cover, hold up a strong magnet - without touching the glass - until the motion of the system is reversed.

The experiment setup must not be subjected to any shocks or temperature variations which may cause convection in the torsion balance housing while recording measured values (see Fig. 6)

Check the zero-point stability before starting the measurement. If necessary, observe and note the zero-point fluctuations for at least 10 minutes so that you can determine an average value for x_0 .

At time $t = 0$ rotate the carrier with the lead balls quickly from one extreme position to the other. However, in carrying out this movement, make sure that fingers or lead balls do not come into contact with the housing. Start the stop-clock immediately after shifting the carrier.

Read off and note the position of the light pointer on the scale every 30 s (Fig. 6) for at least 3 periods when using the end displacement method and for 1 period when using the acceleration

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24-hour log of the position of the torsion pendulum: position changes are essentially due to (day)time dependent variations in temperature (particularly noticeable for exposure to direct sunlight)

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Oscillations of the gravitation torsion balance around the final equilibrium position x_{00}

method.

Measuring example:

Equipment constants:

- Mass of the large lead ball: $m_1=1.5 \text{ kg}$
- Distance between ball center and axis of rotation: $d = 0.05 \text{ m}$
- Distance between center point of the large ball (when touching housing) and the small ball (in equilibrium position): $b = 0.047 \text{ mm}$

Distance between mirror and torsion balance: $L = 4.425 \text{ m}$

Oscillation period T (from Fig. 7):

$$T = (2790 - 315)/4 \text{ s} = 618.8 \text{ s}.$$

Equilibrium position at start: $x_0 = 47 \text{ cm}$

Equilibrium position at end x_{00} (determined from three sequential extremes) e.g.

$$x_{00} = \frac{(x_1+x_3)/2 + x_2}{2} = \frac{x_1}{4} + \frac{x_2}{2} + \frac{x_3}{4} = 62.3 \text{ cm}$$

Difference S between light pointer positions for the initial and final pendulum equilibrium states:

$$S = x_{00} - x_0 = 62.3 \text{ cm} - 47 \text{ cm} = 15.3 \text{ cm}$$

Substitution in equation (1) for the end displacement method gives the gravitation constant f as

$$f = \frac{\pi^2 \cdot (0.047 \text{ m})^2 \cdot 0.05 \text{ m} \cdot 0.153 \text{ m}}{1.5 \text{ kg} \cdot (618.8 \text{ s})^2 \cdot 4.425 \text{ m}}$$

$$= 6.56 \cdot 10^{-11} \text{ m}^3 \text{ kg}^{-1} \text{ s}^2.$$

3.5 Zero-point adjustment

The objective of zero-point adjustment is to position the torsion pendulum so that the dumbbell-shaped pendulum body ⑤ (see Fig. 2) is aligned parallel to the glass plates of the housing in the equilibrium position.

An adjustment is necessary:

possibly before first using the system if the manufacturer's adjustment has been modified, e.g. as a result of transportation, after improper handling (through uncontrolled rotation of torsion head) ③),

after a new torsion band has been inserted as described in section 4.

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Preparation for zero-point adjustment:

Assemble the torsion balance with light pointer arrangement on the lab bench e.g. as described in section 3.3 (or in accordance with the corresponding information in the Instruction Sheet 332 11), but without the large lead balls.

Unlock the torsion pendulum by loosening the screws ⑦ (see Fig. 2).

Using the levelling screws on the stand base (or on the optical bench), align the torsion balance so that the torsion pendulum hangs vertically (check: pin at the end of the pendulum hangs in the axis of the hole).

After these preparatory steps have been carried out, allow the arrangement to hang disturbance-free for around 1 day and then carry out the zero-point adjustment.

Zero-point adjustment:

Lock the torsion pendulum with screw ⑦ and set up the torsion balance at the experiment location (e.g. wall mounting as shown in Fig. 3); then unlock the pendulum to set it in motion.

Note the end deflections of the light pointer on both sides of the scale.

Observe the equilibrium position to which the system tends; to shorten the observation time damp the oscillations diamagnetically: each time one of the pendulum lead balls approaches the glass cover, use a strong magnet to reverse the motion of the system - without touching the glass.

Should the light pointer not show a tendency towards an equilibrium position roughly in the center between the two noted direction reversal points, unscrew the knurled screws (3.1) and turn the torsion head ③ through a small angle in the direction of the "desired zero-point".

Once again observe the two end deflections and the zero-point tendency of the light pointer.

Continue this zero-point adjustment until the "desired zero-point" in the middle between the two direction reversal points of the light pointer has been reached.

4 Replacing the torsion band

A defective torsion band can be replaced by a new band which is supplied in a ready-for-assembly state as spare part (683 21). The replacement is carried out as follows:

Preparation

- Assemble the torsion balance (with locked pendulum body) on a stand base which can be levelled (e.g. 300 01/06) in such a way that the back side is at the front; to ensure safe handling, support the protective tube with replacement torsion band (683 21) in a clamp (Fig. 7).
- Remove screw (a₁); subsequently loosen the other screws fixing the glass plate until it can be removed (Fig. 8.2).
- Remove glass plate (a₂) and cardboard seal (a₃) and pull out the front section of the cardboard slide (a₄).

Removing the defective torsion band

- Loosen grub screw (b₁) while holding the band holder (b) tightly. Remove the band holder (Fig. 8.3).
- If the torsion band is broken, remove end piece (c) (Fig. 8.4).

Inserting the new torsion band

- Loosen screw x₁ and remove screw x₂ while holding band holder (B) tightly.
- Very carefully remove the arrangement from the protective cover and "thread" it into the torsion balance until the end piece (C) hangs 1 mm to 2 mm above the hole (c₂). Make sure the very sensitive torsion band does not come into contact with the walls. Subsequently fix pendulum holder (B) by tightening screw ((b₁) (Fig. 8.6).
- Wait until the torsion band has settled and, if necessary, use the levelling screws on the stand base to align the torsion band so that the end piece (C) hangs vertically above the hole.

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- Loosen screw (b₁). Hold the band holder (B) tightly and lower it a few mm - without twisting it - until the machined part of the end piece (C) has been completely guided into the hole (c₁); do not release the band holder until it has been carefully fixed in place by tightening the grub screw (b₁) (Fig. 8.7); finally fix the end piece (C) with screw (c₁).
- Insert the front section of the slide (a₄). After correctly reinserting the cardboard disc (a₃) and the glass plate (a₂), close the housing with screw (a₁).
- Insert cardboard slide (a₄), and then reclose the housing with screw (a₁) after correctly inserting sealing disc (a₃) and glass plate (a₂).

After assembly has been carried out, unlock the torsion pendulum and check if the pendulum body can oscillate freely; should it lie on the locking springs despite unlocking, slightly raise the band holder (B), with screw (b₁) loosened, and again carefully fix in place.

The torsion pendulum must hang undisturbed for at least 12 hours before the required zero-point adjustment as described in section 3.5 can be carried out.


Fig. 8.1


Fig. 8.3


Fig. 8.5


Fig. 8.2


Fig. 8.4


Fig. 8.7



Fig. 8.6

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To protect the sensitive torsion band during transport, the oscillation system of the gravitation torsion balance has been lowered. Before you can use the balance, you must first carry out the following steps, in addition to the instructions in the instruction sheet.

First, mount the balance without the large lead balls vertically in a stand base.

Completely unscrew the two locking screws ⑦ for the pendulum body (Fig. 2) until the clamps are loosened.

Loosen the grub screw (3.2), Fig. 2, without turning the knurled screw.

Hold onto the band holder (b), Fig. 8.3, so that it does not twist.

Carefully withdraw the band holder (b), Fig. 8.3, and observe the motion of the dumbbell-shaped pendulum body ⑤ (Fig. 2). Lift the pendulum body until the two small lead balls are in the middle of the round window. If the band holder is withdrawn so far that the pendulum body makes contact, the torsion band will break.

Try not to twist the band holder, as this requires readjustment of the zero point (see section 3.5 of the instruction sheet).

Retighten the grub screw (3.2), Fig. 2.

If necessary, readjust the zero point.

Refer to the instruction sheet for all other steps required to set up the torsion balance.

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