
Answer #128

The answer is (a): the reading will indicate a greater voltage or a greater charge, as seen in an mpeg video by clicking your mouse on the photograph below.



When a positively charged rod is brought close to the terminal, it pulls additional electrons out of the electroscope toward the terminal, rendering a greater positive charge on the electroscope, as seen.

Archive 7

Question of the Week

Outreach Index Page

Lecture-Demonstration Home Page



For questions and comments regarding the *Question of the Week* contact <u>Dr. Richard E. Berg</u> by e-mail or using phone number or regular mail address given on the <u>Lecture-Demonstration Home Page</u>.